

# **ImageMedix User Manual**

## **Introduction**

**Welcome to ImageMedix, an AI-powered medical image analysis platform designed to help healthcare professionals analyze medical scans, specifically chest X-rays and brain MRI scans. This user manual provides detailed instructions for using the ImageMedix application.**

## **Table of Contents**

- 1. Landing Page**
- 2. Authentication (Sign-up/Sign-in)**
- 3. Home/Dashboard**
- 4. Upload and Scan**
- 5. Viewing Analysis Results**
- 6. History Page**
- 7. Settings**
- 8. About Page**

## **1. Landing Page**

**When you first access ImageMedix, you'll be greeted by the landing page which includes:**

- Navigation header with links to Features, How It Works, Get Started, and Contact**
- Hero section with a brief description of the platform and call-to-action buttons**
- Features section highlighting capabilities like Rapid Diagnostics, High Accuracy, and Unified Platform**
- How It Works section showing the four-step process: Upload, Analysis, Results, and Decision Making**
- Vision section explaining the technological objectives and future applications**
- Call-to-action section to get started or contact for more information**

**From the landing page, you can:**

- Click "Sign In" in the top right to access your account**
- Click "Get Started" to begin using the platform (redirects to sign-in if not authenticated)**

## **2. Authentication**

## **Sign In**

- 1. Click "Sign In" from the landing page**
- 2. You'll see a sign-in form with:**
  - **Option to sign in with Google (single click)**
  - **Email field for email-based login**
  - **Password field (for email login)**
  - **"Sign up" link if you don't have an account**
- 3. For Google sign-in:**
  - **Click the "Continue with Google" button**
  - **Select your Google account in the popup window**
  - **Authorize the application if prompted**
- 4. For email sign-in:**
  - **Enter your email and password**
  - **Click the sign-in button**
- 5. The system will briefly show a "Completing login..." screen**
- 6. If authentication is successful, you'll be redirected to the Home/Dashboard page**

## **Sign Up**

- 1. From the Sign In page, click "Sign up" at the bottom**
- 2. Choose between:**
  - **Sign up with Google (single click)**
  - **Create account with email**
- 3. For Google sign-up:**
  - **Click the "Continue with Google" button**
  - **Select your Google account in the popup window**
  - **Authorize the application if prompted**
- 4. For email sign-up:**
  - **Fill in your name, email, and password**
  - **Click "Sign up"**
- 5. You'll see a "Processing login..." screen briefly**
- 6. After successful registration, you'll be redirected to the Home page**

## **3. Home/Dashboard**

**After successful authentication, you'll be directed to the Home page which includes:**

- **Sidebar navigation with links to Dashboard, Upload Scans, History, Settings, and About**
- **Welcome message with your name**
- **Three main action cards:**
  - 1. Upload New Scan: Quick access to upload medical scans**
  - 2. Recent Scans: View your recently uploaded scans**

### **3. Settings: Manage your account settings**

**The sidebar remains consistent across all pages for easy navigation. It also displays your user profile at the bottom with a logout option.**

## **4. Upload and Scan**

- 1. From any page, click "Upload Scans" in the sidebar or select "Upload New Scan" from the Home page**
- 2. On the Upload page:**
  - a. Select Scan Type:**
    - **Choose between "Brain Scan" or "Chest X-Ray" using the toggle buttons**
  - b. Enter Patient Information:**
    - **Patient ID (required)**
    - **Patient Name (required)**
  - c. Upload Files:**
    - **Drag and drop files into the upload area or click to browse**
    - **Supported formats: PNG, JPG, JPEG, GIF, DICOM (max 10MB)**
    - **Uploaded files appear in a list below with status indicators**
  - d. Start Analysis:**
    - **Click "Upload & Analyze X-rays" or "Upload Scans" button**
    - **The system will begin processing your files, showing status updates:**
      - **Pending**
      - **Uploading**
      - **Analyzing**
      - **Success/Error**
  - e. Processing:**
    - **Files are analyzed one by one**
    - **Status is updated in real-time**
    - **For each completed analysis, you'll see:**
      - **Condition detected (Normal/Pneumonia for chest X-rays, Normal/Tumor for brain scans)**
      - **Confidence percentage**
      - **"View Report" button**
- 3. Once all files have been analyzed, a success message appears, and you can either:**

- **View detailed results by clicking "View Report" for any individual scan**
- **Click "View All Results" to go to the History page**

## **5. Viewing Analysis Results**

**After a scan is analyzed, you can view detailed results in two ways:**

### **Direct Viewing (from Upload page)**

- 1. Click "View Report" on any successfully analyzed scan**
- 2. A modal window appears with detailed information:**
  - a. Patient Information section:**
    - **Patient name**
    - **Patient ID**
    - **Scan type**
    - **File name**
  - b. Diagnostic Summary section:**
    - **Finding (Normal/Pneumonia or Normal/Brain Tumor)**
    - **Confidence percentage with visual indicator**
    - **Analysis date and time**
  - c. Detailed Report section:**
    - **Technical analysis text**
    - **Option to download the report as a text file**
  - d. AI Analysis section:**
    - **Explanation of findings from the AI system**
  - e. Actions:**
    - **Download Report button**
    - **Close button**

### **From History Page**

- **Results can also be viewed from the History page (detailed in section 6)**

## **6. History Page**

**The History page displays all previously uploaded scans and their results:**

- 1. Access the History page by clicking "History" in the sidebar or "Recent Scans" on the Home page**
- 2. The page displays a list of all your uploaded scans with:**
  - **Patient ID**
  - **Upload date and time**
  - **Scan type**
  - **Status (completed, processing)**
  - **Diagnosis result (if completed)**
  - **Confidence percentage (if completed)**
- 3. For each scan, you have two options:**
  - **View Details: Opens a detailed view modal similar to the one in the Upload page**
  - **Delete: Removes the scan from your history (with confirmation)**
- 4. When viewing scan details:**
  - **See all diagnostic information**
  - **Download the report as a text file**
  - **View the original scan image**
  - **Close to return to the History list**
- 5. The page also includes:**
  - **"Upload New Scan" button to quickly add more scans**
  - **Refresh button to update the list**

## **7. Settings**

**The Settings page allows you to customize your account:**

- 1. Access Settings by clicking "Settings" in the sidebar or the Settings card on the Home page**
- 2. The page is divided into three sections:**
  - a. Profile Settings:**
    - **Full Name (editable)**
    - **Email (non-editable, linked to your account)**
  - b. Preferences:**
    - **Email Notifications toggle (receive notifications about scan analysis results)**
    - **Dark Mode toggle (use dark theme for better visibility)**
  - c. Security:**
    - **Change Password link**
- 3. After making changes, click "Save Changes" to update your settings**

### **Change Password**

1. Click "Change Password" in the Security section
2. Fill out the form:
  - Current Password
  - New Password
  - Confirm New Password
3. Click "Update Password" to save changes
4. You'll be redirected back to the Settings page after successful update

## 8. About Page

The About page provides information about ImageMedix, including:

- Company information
- Technology details
- Contact information
- Team members

## Additional Features

### Logout

- Click "Logout" at the bottom of the sidebar to sign out
- You'll be redirected to the Sign In page

### Viewing Downloaded Reports

- Reports are downloaded as text files (.txt)
- You can open them with any text editor
- They contain a detailed analysis of the scan with findings and recommendations

## Troubleshooting

### Common Issues

1. Failed uploads:
  - Check your internet connection
  - Ensure file is under 10MB and in a supported format
  - Try refreshing the page and uploading again
2. Analysis errors:
  - The system provides mock data if the actual analysis fails
  - Check that your image is clear and properly oriented
  - Try a different image if issues persist
3. Session timeout:
  - If inactive for too long, you may be logged out

- **Simply sign in again to continue**