

ImageMedix Deployment Manual

Introduction

This deployment manual provides comprehensive instructions for setting up the ImageMedix application infrastructure from the repository at <https://github.com/htmw/2025S-The-Minions/wiki>. The application consists of three main components:

1. Node.js Backend (deployed on Heroku)
2. Machine Learning Backend (deployed separately)
3. Next.js Frontend (deployed on Vercel)

Table of Contents

1. Prerequisites
2. Backend Deployment (Heroku)
3. Machine Learning Backend Setup
4. Frontend Deployment (Vercel)
5. Environment Configuration
6. Testing the Deployment
7. Maintenance and Monitoring

1. Prerequisites

Before beginning deployment, ensure you have:

- Heroku account
- Vercel account
- Git installed locally
- Node.js (v16 or later) installed locally
- npm or yarn package manager
- Access to the ImageMedix repository (<https://github.com/htmw/2025S-The-Minions/wiki>)

2. Backend Deployment (Heroku)

2.1 Prepare the Backend for Deployment

1. Clone the repository:
2. `git clone https://github.com/htmw/2025S-The-Minions.git`
3. `cd 2025S-The-Minions`
4. Create a new package.json file in the root of the server directory if it doesn't exist:
5.

```
{
```
6.

```
  "name": "imagemedix-backend",
```
7.

```
  "version": "1.0.0",
```

```

8.   "description": "ImageMedix Backend API",
9.   "main": "index.js",
10.  "scripts": {
11.    "start": "node index.js",
12.    "dev": "nodemon index.js"
13.  },
14.  "dependencies": {
15.    "express": "^4.18.2",
16.    "cors": "^2.8.5",
17.    "dotenv": "^16.0.3",
18.    "mongoose": "^7.0.0",
19.    "multer": "^1.4.5-lts.1",
20.    "axios": "^1.3.4",
21.    "jsonwebtoken": "^9.0.0",
22.    "bcrypt": "^5.1.0"
23.  },
24.  "devDependencies": {
25.    "nodemon": "^2.0.20"
26.  }
27.}
28.   Create a Procfile in the server directory:
29.web: npm start
30.   Create a .env file in the server directory for local
    development:
31.NODE_ENV=development
32.PORT=8080
33.MONGODB_URI=mongodb://localhost:27017/imagemedix
34.JWT_SECRET=your_local_jwt_secret
35.ML_API_URL=http://localhost:5000
36.CORS_ORIGIN=http://localhost:3000

```

2.2 Deploy to Heroku

```

1. Install the Heroku CLI and log in:
2. npm install -g heroku
3. heroku login
4. Create a new Heroku app:
5. cd server
6. heroku create imagemedix-backend
7. Add the MongoDB add-on for database storage:
8. heroku addons:create mongodb:hobby-dev
9. Configure environment variables:
10.heroku config:set NODE_ENV=production
11.heroku config:set JWT_SECRET=your_production_jwt_secret
12.heroku config:set ML_API_URL=https://your-ml-api-url.herokuapp.com
13.heroku config:set CORS_ORIGIN=https://imagemedix.vercel.app
14.   Deploy the backend to Heroku:
15.git subtree push --prefix server heroku main

```

If the above command fails, you can try:

```
git push heroku `git subtree split --prefix server main`:main --force
```

```

16.   Ensure at least one instance is running:
17.heroku ps:scale web=1

```

18. Check the logs to verify the deployment:
19. `heroku logs --tail`

3. Machine Learning Backend Setup

3.1 Prepare the ML Service

1. Navigate to the ML service directory:
2. `cd ml-service`
3. Create a requirements.txt file if it doesn't exist:
4. `flask==2.2.3`
5. `gunicorn==20.1.0`
6. `numpy==1.24.2`
7. `tensorflow==2.11.0`
8. `pillow==9.4.0`
9. `scikit-learn==1.2.2`
10. `opencv-python-headless==4.7.0.72`
11. Ensure the Dockerfile exists and contains:
12. `FROM python:3.10-slim`
- 13.
14. `WORKDIR /app`
- 15.
16. `COPY requirements.txt .`
17. `RUN pip install --no-cache-dir -r requirements.txt`
- 18.
19. `COPY . .`
- 20.
21. `EXPOSE 5000`
- 22.
23. `CMD ["gunicorn", "--bind", "0.0.0.0:5000", "app:app"]`

3.2 Deploy ML Service to Heroku

1. Create a new Heroku app for the ML service:
2. `heroku create imagemedix-ml-service`
3. Log in to Heroku Container Registry:
4. `heroku container:login`
5. Build and push the Docker container to Heroku:
6. `heroku container:push web --app imagemedix-ml-service`
7. Release the container:
8. `heroku container:release web --app imagemedix-ml-service`
9. Configure environment variables:
10. `heroku config:set MODEL_PATH=/app/models --app imagemedix-ml-service`
11. `heroku config:set ALLOWED_ORIGINS=https://imagemedix-backend.herokuapp.com --app imagemedix-ml-service`
12. Check the logs to verify deployment:
13. `heroku logs --tail --app imagemedix-ml-service`

4. Frontend Deployment (Vercel)

4.1 Prepare the Frontend

1. Navigate to the app directory:
2. `cd app`
3. Install dependencies:
4. `npm install`
5. Create a `.env.local` file:
6. `NEXT_PUBLIC_API_URL=https://imagemedix-backend.herokuapp.com/api`
7. `NEXT_PUBLIC_ML_API_URL=https://imagemedix-ml-service.herokuapp.com`
8. `CLERK_PUBLISHABLE_KEY=your_clerk_publishable_key`
9. `CLERK_SECRET_KEY=your_clerk_secret_key`
10. Verify the build works locally:
11. `npm run build`
12. `npm run start`

4.2 Deploy to Vercel

1. Install the Vercel CLI:
2. `npm install -g vercel`
3. Log in to Vercel:
4. `vercel login`
5. Deploy to Vercel:
6. `vercel`

Follow the prompts to set up the project:

- Set the root directory to `app`
 - Configure the build settings
 - Add environment variables
7. For production deployment:
 8. `vercel --prod`
 9. Alternatively, connect your GitHub repository to Vercel:
 - Go to <https://vercel.com/new>
 - Import your GitHub repository
 - Configure the project settings (set root directory to `app`)
 - Add environment variables
 - Deploy

5. Environment Configuration

5.1 Backend Environment Variables

Variable Name	Description	Example Value
NODE_ENV	Environment mode	production
PORT	Port the server runs on	8080
MONGODB_URI	MongoDB connection string	mongodb://user:pass@host:port/db

JWT_SECRET	Secret for JWT tokens	your_jwt_secret_key
ML_API_URL	URL for the ML service	https://imagedix-ml-service.herokuapp.com
CORS_ORIGIN	Allowed origins for CORS	https://imagedix.vercel.app

5.2 ML Service Environment Variables

Variable Name	Description	Example Value
MODEL_PATH	Path to model files	/app/models
ALLOWED_ORIGINS	Allowed origins for CORS	https://imagedix-backend.herokuapp.com
LOG_LEVEL	Logging level	info

5.3 Frontend Environment Variables

Variable Name	Description	Example Value
NEXT_PUBLIC_API_URL	Backend API URL	https://imagedix-backend.herokuapp.com/api
NEXT_PUBLIC_ML_API_URL	ML service URL	https://imagedix-ml-service.herokuapp.com
CLERK_PUBLISHABLE_KEY	Clerk authentication key	pk_test_XXXXX
CLERK_SECRET_KEY	Clerk authentication secret	sk_test_XXXXX

6. Testing the Deployment

6.1 Backend Testing

1. Test the health endpoint:
2. `curl https://imagedix-backend.herokuapp.com/api/health`
3. Test user registration:
4. `curl -X POST https://imagedix-backend.herokuapp.com/api/auth/register \`
5. `-H "Content-Type: application/json" \`
6. `-d '{"name":"Test User","email":"test@example.com","password":"password123"}'`
7. Test user login:
8. `curl -X POST https://imagedix-backend.herokuapp.com/api/auth/login \`
9. `-H "Content-Type: application/json" \`
10. `-d '{"email":"test@example.com","password":"password123"}'`

6.2 ML Service Testing

1. Test the health endpoint:

2. `curl https://imagemedix-ml-service.herokuapp.com/health`
3. Test the prediction endpoint with a sample chest X-ray:
4. `curl -X POST https://imagemedix-ml-service.herokuapp.com/api/ml/analyze-chest \`
5. `-F "image=@sample-xray.jpg" \`
6. `-H "Content-Type: multipart/form-data"`

6.3 Frontend Testing

1. Open the deployed frontend in a browser:
2. `https://imagemedix.vercel.app`
3. Test all main flows:
 - Registration and login
 - Uploading medical scans
 - Viewing scan analysis results
 - Accessing scan history
 - Updating settings

7. Maintenance and Monitoring

7.1 Backend Monitoring

1. View Heroku logs:
2. `heroku logs --tail --app imagemedix-backend`
3. Set up Heroku metrics:
4. `heroku addons:create librato:development --app imagemedix-backend`
5. Set up error tracking with Sentry:
6. `heroku addons:create sentry:f1 --app imagemedix-backend`

7.2 ML Service Monitoring

1. View ML service logs:
2. `heroku logs --tail --app imagemedix-ml-service`
3. Monitor ML service performance:
4. `heroku addons:create newrelic:wayne --app imagemedix-ml-service`

7.3 Frontend Monitoring

1. View Vercel deployment logs from the Vercel dashboard
2. Set up analytics:
 - Add Vercel Analytics by enabling it in the Vercel dashboard
 - Implement custom analytics using a service like Plausible or Umami

7.4 Scaling Considerations

1. Backend scaling on Heroku:
2. `heroku ps:scale web=2 --app imagemedix-backend`
3. ML service scaling:

4. `heroku ps:scale web=2 --app imagemedix-ml-service`
5. Database scaling:
 - Upgrade MongoDB plan as needed
 - Consider adding read replicas for high traffic

7.5 Backup Procedures

1. Set up automatic MongoDB backups:
2. `heroku addons:create mongolab:backup-daily --app imagemedix-backend`
3. Implement regular model backups:
 - Store ML model files in a separate storage service (AWS S3)
 - Version control your models
 - Document model versions and changes

Conclusion

Following this deployment manual will result in a fully functional ImageMedix application with:

- Node.js backend hosted on Heroku
- Custom ML model service running on a separate Heroku instance
- Next.js frontend deployed on Vercel

Once deployed, the application will be accessible at your custom domain (e.g., `imagemedix.vercel.app`), connecting to your backend (e.g., `imagemedix-backend.herokuapp.com/api`) and utilizing your ML service (e.g., `imagemedix-ml-service.herokuapp.com`) for image analysis.

After completing the deployment, be sure to test the connections between all components thoroughly before making the application available to users.